




MAUDEMetrics: Automated extraction and comprehensive analysis of FDA MAUDE device safety data

Mohamed Marouf ¹

¹ Faculty of Medicine, Mansoura University, Mansoura, Egypt

DOI: [10.xxxxxx/draft](https://doi.org/10.xxxxxx/draft)

Software

- [Review](#) 
- [Repository](#) 
- [Archive](#) 

Editor: 

Submitted: 07 September 2025

Published: unpublished

License

Authors of papers retain copyright and release the work under a Creative Commons Attribution 4.0 International License ([CC BY 4.0](https://creativecommons.org/licenses/by/4.0/)).

Summary

Monitoring medical device safety is critical for patient safety and regulatory compliance. The FDA's Manufacturer and User Facility Device Experience (MAUDE) database contains over 22 million adverse event reports and remains a primary source for postmarket surveillance ([U.S. Food and Drug Administration, 2025](#)).

The FDA receives several hundred thousand new medical device reports annually, which are accessible via the openFDA Device Event API ([openFDA, 2025](#)).

However, practical use of this database is constrained by the API's 1,000-record pagination limit, the FDA website's 500-record export restriction (with missing demographic fields), and the need for significant technical expertise to aggregate and analyze reports.

These barriers make comprehensive device safety analysis difficult, time-consuming, and error-prone.

MAUDEMetrics is a Python-based web application that democratizes access to MAUDE data by automating extraction, aggregation, and analysis of adverse event reports via the openFDA API.

It provides an intuitive web interface, requires no programming expertise, and supports complex queries across brand names, product codes, manufacturers, and date ranges.

The tool automatically handles API pagination, validates and standardizes fields, and produces professional multi-sheet Excel exports with raw data, cleaned datasets, and statistical summaries.

By reducing analysis time from hours to minutes, MAUDEMetrics enhances reproducibility, efficiency, and accessibility of MAUDE data for clinicians, researchers, quality assurance teams, and regulatory professionals.

Statement of Need

Medical device safety analysis is essential for patient safety and regulatory compliance, yet current MAUDE access methods impose significant barriers:

- **API Limitations:** The openFDA API restricts results to 1,000 records per query, requiring complex iteration ([openFDA, 2025](#)).
- **Web Interface Limitations:** The MAUDE portal only allows 500-record exports and omits essential fields such as demographics ([U.S. Food and Drug Administration, 2025](#)).
- **Manual Burden:** Researchers must perform multiple searches, downloads, and merges to assemble datasets.

40 ▪ **Reproducibility Challenges:** Manual extraction is difficult to replicate and error-prone.

41 **MAUDEMetrics** addresses these gaps by providing:

42 ▪ An accessible, programming-free web interface

43 ▪ Intelligent pagination to assemble complete datasets

44 ▪ Multi-parameter search across brand names, product codes, manufacturers, product
45 classes, and date ranges

46 ▪ Automated cleaning and validation of >100 FDA data fields

47 ▪ Multi-sheet professional Excel exports with summaries

48 ▪ Built-in visualization dashboard for event types and device brands

49 ▪ Docker containerization for reproducible workflows

50 The software benefits clinical researchers, hospital QA teams, regulatory professionals, and
51 healthcare institutions conducting device safety audits.

52 **Implementation**

53 MAUDEMetrics is implemented as a **Flask** web application ([Flask Development Team, 2023](#))
54 that integrates directly with the openFDA Device Event API.

55 The `fetch_all_API_data()` function automates API pagination and error handling, while
56 **Pandas** ([McKinney, 2010](#)) and **NumPy** ([Harris et al., 2020](#)) perform data processing.

57 Data are cached locally in **SQLite** ([SQLite Development Team, 2023](#)) for performance and
58 persistence.

59 The export engine provides two complementary options for researchers:

60 ▪ **Raw Data Export:** reports are saved as CSV files that preserve the original JSON
61 structure from the openFDA API, ensuring no loss of information.

62 ▪ **Processed Excel Export:** using **OpenPyXL** ([Clark, 2023](#)), the system generates multi-
63 sheet Excel workbooks with:

64 – **Events:** cleaned and standardized event data with integrated narratives

65 – **Summary:** aggregated statistics of patient demographics and event types

66 Deployment uses **Docker** ([Docker Inc., 2023](#)) for reproducibility, and the web interface is styled
67 with Bootstrap 5.

68 Testing covers API integration, database operations, export generation, and web routes, with
69 a 100% pass rate.

70 Performance benchmarks show ~5,000 records extracted in under 30 seconds.




Figure 1: MAUDEMetrics Workflow

77 *Figure 1: MAUDEMetrics architecture and workflow. User input via the web interface triggers*
 78 *Flask endpoints that query the openFDA API (with pagination), cache results in SQLite,*
 79 *and process data with Pandas. The system supports two export modes: raw CSV files that*
 80 *preserve the original FDA JSON structure, and processed Excel workbooks (via OpenPyXL)*
 81 *with cleaned datasets and summary statistics. Docker ensures reproducible deployment.*

82 Use Cases and Impact

83 Research:

- 84 - Aggregate device-specific adverse events for systematic reviews
- 85 - Compare event rates across manufacturers
- 86 - Conduct temporal trend analyses

87 Clinical Quality Assurance:

- 88 - Monitor devices used within an institution
- 89 - Detect emerging safety signals
- 90 - Generate reproducible safety reports

91 Regulatory:

- 92 - Support FDA submissions with structured evidence
- 93 - Standardize postmarket surveillance workflows

94 By automating extraction and standardizing outputs, MAUDEMetrics reduces manual effort,
 95 minimizes transcription errors, and enables reproducible device safety investigations.

96 Software Availability

- 97 ▪ **Name:** MAUDEMetrics
- 98
- 99 ▪ **Repository:** <https://github.com/MohamedMaroufMD/MAUDEMetrics>
- 100
- 101 ▪ **DOI:** <https://doi.org/10.5281/zenodo.16691960>
- 102
- 103 ▪ **License:** Apache License 2.0
- 104
- 105 ▪ **Version:** 2.1.0 (2025-09-08)

106 Detailed installation (Docker and manual), usage examples, tests, and sample exports are
 107 provided in the repository.

108 References

- 109 Clark, C. (2023). *OpenPyXL: A python library for excel 2010 xlsx/xlsm files*. <https://openpyxl.readthedocs.io/>
- 110
- 111 Docker Inc. (2023). *Docker: Empowering app development for developers*. <https://www.docker.com/>
- 112
- 113 Flask Development Team. (2023). *Flask: A lightweight WSGI web application framework*.
 114 <https://flask.palletsprojects.com/>
- 115 Harris, C. R., Millman, K. J., Walt, S. J. van der, Gommers, R., Virtanen, P., Cournapeau, D.,
 116 Wieser, E., Taylor, J., Berg, S., Smith, N. J., Kern, R., Picus, M., Hoyer, S., Kerkwijk,
 117 M. H. van, Brett, M., Haldane, A., Río, J. F. del, Wiebe, M., Peterson, P., ... Oliphant,
 118 T. E. (2020). Array programming with NumPy. *Nature*, 585(7825), 357–362. <https://doi.org/10.1038/s41586-020-2649-2>
- 119

- 120 McKinney, W. (2010). Data structures for statistical computing in python. *Proceedings of the*
121 *9th Python in Science Conference*, 51–56. <https://doi.org/10.25080/majora-92bf1922-00a>
- 122 openFDA. (2025). *openFDA device event API*. <https://open.fda.gov/apis/device/event/>
- 123 SQLite Development Team. (2023). *SQLite: Self-contained, high-reliability SQL database*
124 *engine*. <https://www.sqlite.org/>
- 125 U.S. Food and Drug Administration. (2025). *Manufacturer and user facility device experi-*
126 *ence (MAUDE)*. [https://www.fda.gov/medical-devices/postmarket-requirements-devices/](https://www.fda.gov/medical-devices/postmarket-requirements-devices/mandatory-reporting-requirements-manufacturers-importers-and-device-user-facilities)
127 [mandatory-reporting-requirements-manufacturers-importers-and-device-user-facilities](https://www.fda.gov/medical-devices/postmarket-requirements-devices/mandatory-reporting-requirements-manufacturers-importers-and-device-user-facilities)

DRAFT