




MacroQA: An ImageJ Macro for ACR MRI Quality Assurance

Gabriel Branco Vitorino^{1*}, Victor Hugo Celoni Gnatkovski^{1*}, Pedro Henrique Tosta Cayres de Oliveira^{1*}, Mateus Setubal Quiel^{1*}, and Ícaro Agenor Ferreira de Oliveira ^{2,3¶}

¹ Departamento de Física, Faculdade de Filosofia, Ciências e Letras de Ribeirão Preto, Universidade de São Paulo, Brazil ² Fundação de Apoio ao Ensino, Pesquisa e Assistência do Hospital das Clínicas da Faculdade de Medicina de Ribeirão Preto, Brazil ³ Centro das Ciências da Imagem e Física Médica, Hospital das Clínicas da Faculdade de Medicina de Ribeirão Preto, Universidade de São Paulo, Brazil ¶ Corresponding author * These authors contributed equally.

DOI: [10.xxxxxx/draft](https://doi.org/10.xxxxxx/draft)

Software

- [Review](#) 
- [Repository](#) 
- [Archive](#) 

Editor: 

Submitted: 10 October 2025

Published: unpublished

License

Authors of papers retain copyright and release the work under a Creative Commons Attribution 4.0 International License ([CC BY 4.0](https://creativecommons.org/licenses/by/4.0/)).

Summary

MacroQA is an open-source ImageJ/Fiji macro package that implements the American College of Radiology (ACR) quality assurance (QA) tests for MRI phantoms. The project was developed with academic and pedagogical goals in mind, and it aims to simplify and standardize phantom testing. By leveraging Fiji/ImageJ's built-in functionality, MacroQA performs the ACR phantom tests quickly and reproducibly, completing the QA workflow within minutes. As a free and accessible alternative to proprietary software, MacroQA lowers barriers to adoption, promotes reproducibility, and supports collaborative development in the MRI research and clinical communities.

Statement of need

Magnetic resonance imaging (MRI) is indispensable in modern medicine and neuroscience, enabling diverse applications from clinical diagnosis to the study of functional brain connectivity. (Dumoulin et al., 2018; Granziera et al., 2021; MacDonald & Frayne, 2015; Rüber et al., 2018). To uphold image quality and ensure the reproducibility of results—a central concern across the MR community—a robust Quality Assurance (QA) and Quality Control (QC) procedures are essential (Epistatou et al., 2020; Stöcker et al., 2025; Sun et al., 2015; Vogelbacher et al., 2019). The most widely adopted standards is the QA program established by the American College of Radiology (ACR), which relies on a dedicated accreditation phantom. Despite the critical need for standardization, existing software solutions for ACR QA often depend on proprietary platforms. These require expensive licenses and operate within closed-source ecosystems, creating significant financial and accessibility barriers for research, educational, and clinical facilities worldwide. MacroQA directly addresses this fundamental gap. By providing a comprehensive, open-source implementation of ACR QA test suite within Fiji/ImageJ platform, MacroQA eliminates the requirement for commercial software. Implemented in the Jython scripting language, this design promotes transparency, verifiability, and accessibility, offering a cost-effective, shareable solution that supports both clinical best practices and reproducible research.

Installation

1. Ensure that you have [Fiji](#) installed, preferably with Java 8 runtime. *Note: We recommend using the Fiji distribution because it already includes the Jython library.*
2. Clone or download the MacroQA repository from this [GitHub](#) page. *Note: This software is a self-contained ImageJ/Fiji macro and does not require any external dependencies beyond a standard installation of Fiji. It relies solely on the core functions of ImageJ and Jython.*

How to use MacroQA in Fiji?

MacroQA can be used in two main ways, depending on your preference:

Method 1: Run directly via Macro Editor

This method is ideal for quick use or one-off tests.

Steps (Fiji): 1. Open the **StartupMacros** in the *Plugins > Macros* tab. 2. In your file explorer, locate the MacroQA folder. 3. Open the folder and double-click on the macro that you want to run. 4. The macro will open in Fiji's editor - simply press *Run*.

For ImageJ/ImageJ2 users the steps are similar, but ensure that the Jython library is also installed.*

Method 2: Install as a Plugin

Installing MacroQA as a plugin integrates it into Fiji's menu system, making it persistently available across sessions.


Steps (general): 1. Copy the MacroQA folder into a subdirectory of your Fiji plugins folder (for example, `.../Fiji.app/plugins/` or `.../Fiji.app/plugins/Macros/`). 2. Restart Fiji. 3. The macros will now appear in the *Plugins > Macros* menu.

Platform-specific examples: - Windows (typical): `C:\Program Files\Fiji\Fiji.app\plugins\Macros\MacroQA`
- macOS (typical): `/Applications/Fiji.app/plugins/Macros/MacroQA` or `~/Fiji.app/plugins/Macros/MacroQA`
- Linux (typical): `/home/<user>/Fiji.app/plugins/Macros/MacroQA` or `/opt/Fiji.app/plugins/Macros/MacroQA`


Functionality

Usage example


Once installed, MacroQA becomes available under the *Plugins > Macros* menu in Fiji (Figure 1). From there, the user can select any of the available ACR quality control tests, such as Central Frequency, Geometric Accuracy, or Signal-to-Noise Ratio.



68 When a test is launched, MacroQA
69 guides the user through the required steps via dialog boxes and messages. For example,
70 running the Central Frequency test (Figure 2) prompts the user to select the appropriate
71 image series and automatically reports the measured resonance frequency in the Fiji log



72 window. Some tests require user
73 interaction, such as drawing straight lines or selecting regions of interest. In the Geometric
74 Accuracy test (Figures 3), the macro requests that the user load the Localizer image
75 and draw reference lines across the phantom. It then requests that the user load the
76 ACR T1 series, where two different slices are assessed. These inputs are then used to
77 calculate geometric dimensions, which are compared against the ACR acceptance criteria.



This combination of guided prompts and automated calculations ensures that even users with limited prior experience can reliably perform ACR phantom quality control tests in a reproducible manner.

Quality control tests and their acceptance criteria

Required images

For both large and medium phantoms, a minimum of three acquisitions are required: the **Localizer**, **ACR T1 series**, and the **ACR T2 series**. - **Localizer**: a single-slice sagittal spin-echo acquired at the phantom's center. - **ACR T1**: an 11-slice axial T1-weighted (T1w) series. - **ACR T2**: an 11-slice axial T2-weighted (T2w) series acquired with two echo times; the longer echo is used as the T2-weighted image.

Below is a brief summary of the quality control tests supported by MacroQA. Users are encouraged to first review and follow the [ACR MRI Phantom testing guidelines](#) when using MacroQA for the first time. This ensures familiarity with the procedures and acceptance criteria before relying on automated analysis.

Central frequency

Objective: Ensure the scanner operates at the correct resonance frequency. Off-resonance operation reduces signal-to-noise ratio (SNR) and may indicate drift in the static magnetic field.

Frequency: weekly

Acceptance criteria: within 1 ppm per day for superconducting magnets

Image type: ACR T1-weighted (T1w)

Geometric accuracy

Objective: Verify that image scaling reflects the true dimensions of the imaged object.

102 **Frequency:** weekly

103 **Acceptance criteria:** ± 3 mm (large phantom) and ± 2 mm (medium phantom)

104 **Image type:** ACR T1-weighted (T1w)

105 High-contrast spatial resolution

106 **Objective:** Assess the scanner's ability to resolve small objects.

107 **Frequency:** weekly

108 **Acceptance criteria:** visualization of the 1 mm holes

109 **Image type:** ACR T1-weighted (T1w) and T2-weighted (T2w)

110 Slice thickness accuracy

111 **Objective:** Verify that the prescribed slice thickness matches the acquired slice.

112 **Frequency:** annual

113 **Acceptance criteria: ± 0.75 mm**

114 **Image type:** ACR T1-weighted (T1w) and T2-weighted (T2w)

115 Slice position accuracy

Objective: Assess the accuracy of slice positioning using the localizer image as a reference.

117 **Frequency:** annual

118 **Acceptance criteria:** 5 mm in both directions

119 **Image type:** ACR T1-weighted (T1w) and T2-weighted (T2w)

120 **Image intensity uniformity**

121 **Objective:** Measure intensity uniformity over a large water-only region of the phantom near
122 the middle of the imaged volume (typically near the head coil center).

123 **Frequency:** annual

Acceptance criteria: for scanners at 3T: PIU ≤ 80 ; for scanners < 3 T: PIU ≤ 85

125 **Image type:** ACR T1-weighted (T1w) and T2-weighted (T2w)

126 Percent-signal ghosting

127 **Objective:** Quantify ghosting artifacts in ACR images.

128 **Frequency:** annual

129 **Acceptance criteria:** 3%

130 **Image type:** ACR T1-weighted (T1w) and T2-weighted (T2w)

131 Low-contrast object detectability

132 **Objective:** Determine the extent to which low-contrast objects are discernible in the images.

133 **Frequency:** weekly

134 **Acceptance criteria:** for scanners at 3T: 37 spokes (ACR T1 and T2). For scanners between
135 1.5T and <3T: 30 spokes (ACR T1) and 25 spokes (ACR T2).

136 **Image type:** ACR T1-weighted (T1w) and T2-weighted (T2w)

137 **Signal-to-noise ratio (SNR)**

138 **Objective:** Measure the ratio of true signal to background noise. Although SNR is not always
139 explicitly included in the ACR manual, it is a key indicator of image quality.

140 **Frequency:** weekly

141 **Acceptance criteria:** not formally specified by the ACR

142 **Image type:** ACR T1-weighted (T1w) — note: SNR may require additional acquisitions or
143 specific measurement regions

144 **Software description**

145 MacroQA is implemented in Jython (www.jython.org), the Python implementation for the
146 Java platform, and runs within Fiji/ImageJ. Fiji was chosen because it is widely used, free,
147 and cross-platform. Each QC test is implemented as an independent macro, which simplifies
148 development and installation. After installation, the macros appear in Fiji's menu under a
149 dedicated "MacroQA" submenu.

- 150 ■ **Inputs:** DICOM images acquired with the ACR accreditation phantom.
- 151 ■ **Outputs:** numerical results displayed in the Fiji log window and optionally saved to disk.

152 MacroQA is distributed under the GNU General Public License v3.0 (GPL-3.0), which ensures
153 the code remains free to use, modify, and redistribute under the license terms.

154 **Availability**

155 MacroQA is publicly available on [GitHub](https://github.com). The tool is distributed under the GNU General
156 Public License v3.0 (GPL-3.0). Installation and usage instructions are provided in the repository
157 README. We welcome contributions and feedback from the community — please open an
158 issue to report bugs or request features; pull requests are also welcome.

159 **Acknowledgements**

160 GBV, VHCG, PHTCO, and MSQ (students in the Medical Physics program at the University
161 of São Paulo) contributed nearly equally to coding, software development, and manuscript
162 preparation during their internship. IAFO conceived the idea for MacroQA, contributed
163 to coding and writing, and provided supervision and mentorship throughout the students'
164 internship at the MRI service of the Ribeirão Preto Medical School, University of São Paulo.

165 **References**

- 166 Dumoulin, S. O., Fracasso, A., Zwaag, W. van der, Siero, J. C. W., & Petridou, N. (2018).
167 Ultra-high field MRI: Advancing systems neuroscience towards mesoscopic human brain
168 function. *NeuroImage*, 168, 345–357. <https://doi.org/10.1016/j.neuroimage.2017.01.028>
- 169 Epistatou, A. C., Tsalafoutas, I. A., & Delibasis, K. K. (2020). An automated method for
170 quality control in MRI systems: Methods and considerations. *Journal of Imaging*, 6(10),
171 111. <https://doi.org/10.3390/jimaging6100111>
- 172 Granziera, C., Wuerfel, J., Barkhof, F., Calabrese, M., De Stefano, N., Enzinger, C., Evangelou,
173 N., Filippi, M., Geurts, J. J. G., Reich, D. S., Rocca, M. A., Ropele, S., Rovira, À., Sati,
174 P., Toosy, A. T., Vrenken, H., Gandini Wheeler-Kingshott, C. A. M., & Kappos, L. (2021).
175 Quantitative magnetic resonance imaging towards clinical application in multiple sclerosis.
176 *Brain*, 144(5), 1296–1311. <https://doi.org/10.1093/brain/awab029>

- 177 MacDonald, M. E., & Frayne, R. (2015). Cerebrovascular MRI: A review of state-of-the-
178 art approaches, methods and techniques. *NMR in Biomedicine*, 28(7), 767–791. <https://doi.org/10.1002/nbm.3322>
179
- 180 Rüber, T., David, B., & Elger, C. E. (2018). MRI in epilepsy: Clinical standard and
181 evolution. *Current Opinion in Neurology*, 31(2), 223. [https://doi.org/10.1097/WCO.](https://doi.org/10.1097/WCO.0000000000000539)
182 [0000000000000539](https://doi.org/10.1097/WCO.0000000000000539)
- 183 Stöcker, T., Keenan, K. E., Knoll, F., Priovoulos, N., Uecker, M., & Zaitsev, M. (2025).
184 Reproducibility and quality assurance in MRI. *Magnetic Resonance Materials in Physics,*
185 *Biology and Medicine*, 38(3), 347–352. <https://doi.org/10.1007/s10334-025-01271-1>
- 186 Sun, J., Barnes, M., Dowling, J., Menk, F., Stanwell, P., & Greer, P. B. (2015). An open
187 source automatic quality assurance (OSAQA) tool for the ACR MRI phantom. *Australasian*
188 *Physical & Engineering Sciences in Medicine*, 38(1), 39–46. [https://doi.org/10.1007/](https://doi.org/10.1007/s13246-014-0311-8)
189 [s13246-014-0311-8](https://doi.org/10.1007/s13246-014-0311-8)
- 190 Vogelbacher, C., Bopp, M. H. A., Schuster, V., Herholz, P., Jansen, A., & Sommer, J. (2019).
191 LAB-QA2GO: A free, easy-to-use toolbox for the quality assessment of magnetic resonance
192 imaging data. *Frontiers in Neuroscience*, 13. <https://doi.org/10.3389/fnins.2019.00688>

DRAFT